

Side Shelf

The side shelf can be mounted to either side of the upper side rails of the ScienceDesk, providing a shelf width of 600mm (23.6").


The top surface of the shelf is coated with electrostatic-dissipative material to minimize electrostatic discharge. The shelf has a 75kg (165lbs) load capacity.



Side Shelf: Stocked for Immediate Delivery

ITEM#	\$	£	€	RMB	DESK LENGTH	SHELF DIMENSION
PSX352	\$ 389.30	£245.30	€ 362.00	¥3,717.80	600mm	605mm x 600mm (23.8" x 23.6")
PSX350	\$ 396.10	£249.50	€ 368.40	¥3,782.80	750mm	755mm x 600mm (29.7" x 23.6")
PSX351	\$ 397.80	£250.60	€ 370.00	¥3,799.00	900mm	905mm x 600mm (35.6" x 23.6")

PC Moveable Shelf



The PC shelf can be mounted to the upper or lower rails of the ScienceDesk. It provides a 450mm x 260mm (17.7" x 10.2") shelf size.

The top surface is coated with electrostatic-dissipative material to minimize electrostatic discharge, and it has a 40kg (88lbs) load capacity.

The shelf is ergonomically positioned on the frame of the ScienceDesk™, which is ideally suited for supporting a PC.

PC/Moveable Shelf: Stocked for Immediate Delivery

ITEM#	\$	£	€	RMB	SHELF DIMENSION
PSX180	\$ 255.00	£ 160.70	€ 237.20	¥ 2,435.30	450mm x 260mm (17.7" x 10.2")

Please contact your local sales office for a quotation including shipping.

ScienceDesk™ Keyboard Holders

The frame-mounted Keyboard Holder is suitable for most keyboards and mouse pads. It provides 360° horizontal positioning.

Keyboard Holders: Stocked for Immediate Delivery

ITEM#	\$	£	€	RMB	DESCRIPTION
PSX132	\$ 480.00	£302.40	€446.40	¥ 4584.00	Frame Mounted Keyboard Support




TOOLS OF THE TRADE

Air Compressors

For laboratories not equipped with a compressed air supply, Thorlabs provides a low-noise air compressor that is self-contained, lightweight, and compact. The compressor will not interfere with sensitive experiments because of its low acoustic noise design (35dB at 1m). Each system includes its own air regulator, filter, and receiver tank with pressure gauge, safety valve, drain, and line pressure gauge.

See Page 68 for More Details.



Pressure Control Panel

The Pressure Control Panel allows the passive isolation system to be supplied from a constant air supply, such as a gas bottle or compressor.

See Page 68 For More Details

