


Tables & Breadboards

ScienceDesk™ Adjustable Wrist Rest

The ergonomic wrist rest, designed to provide the optimum height and angle for adjusting nanopositioning stages, enables operators to perform delicate instrument adjustments with higher accuracy and lower fatigue. The adjustable wrist rest mounts to the upper front bar of the ScienceDesk™ base. The top surface is coated with electrostatic-dissipative material to minimize electrostatic discharge. It comes complete with integral mounting studs and locking bolts.



Adjustable Wrist Rest: Stocked for Immediate Delivery

ITEM#	\$	£	€	RMB	FOR BREADBOARDS	DIMENSION
PSX100	\$ 399.50	£ 251.70	€ 371,50	¥ 3,815.20	60mm (2.4") Thick or 55mm (2.2") Thick	800 x 60mm (31.5" x 2.4")
PSX101	\$ 399.50	£ 251.70	€ 371,50	¥ 3,815.20	110mm (4.3") Thick	800 x 60mm (31.5" x 2.4")

Flat Screen Support Bracket



The flat screen support bracket is equipped with a multidirectional pivot clamp, which supports most types of flat screens.


The mounting post (included) is 600mm (23.6") high.



Flat Screen Support Bracket: Stocked for Immediate Delivery

ITEM#	\$	£	€	RMB	DESCRIPTION
PSX120	\$ 323.00	£ 203.50	€ 300,40	¥ 3,084.70	Flat Screen Support Bracket

Post Assembly




These post assemblies are suitable for mounting auxiliary equipment. Three models are offered: the first with a 42mm (1.7") outside diameter and a length of 600mm (23.6"), the second with a 25mm (1.0") outer diameter and a length of 200mm (7.9"), and the third 25mm (1.0") outer diameter and a length of 600mm (23.6").

Post Assemblies: Call for Lead Time

ITEM#	\$	£	€	RMB	SIZE
PSX150	\$ 93.50	£ 58.90	€ 87,00	¥ 892.90	Ø25 x 200mm
PSX151	\$ 110.50	£ 69.60	€ 102,80	¥ 1,055.30	Ø25 x 600mm
PSX160	\$ 127.50	£ 80.30	€ 118,60	¥ 1,217.60	Ø42 x 600mm

Earth Grounding Kit

The Earth Grounding Kit contains a variety of straps, wires, cables, and fittings to enable the user to connect the breadboard, frame, legs, and shelves to ensure proper earth grounding of the entire system.



PSX105



Earth Grounding Kit: Stocked for Immediate Delivery

ITEM#	\$	£	€	RMB	DESCRIPTION
PSX105	\$ 168.30	£ 106.00	€ 156,50	¥ 1,607.30	Frame Earth Grounding Kit

Earthing Kits should be ordered at the same time as the ScienceDesk.