

For current pricing, please see our website.

CHAPTERS

Tables/
Breadboards

Mechanics

Optomechanic
Devices

Kits

Lab Supplies

SECTIONS

Breadboards

Breadboard
Supports

Breadboard
Accessories

ScienceDesk

ScienceDesk
Accessories

Optical Tables

Table Supports

Optical Table
Accessories


Free-Standing Instrumentation Shelf

The Thorlabs free-standing instrumentation shelf provides storage space or further mounting surfaces for instruments, power supplies, and other accessories. The system is available in two sizes that span 6.0' (2.0 m) and 8.0' (2.5 m). Simply match the shelf length to the length or width of your optical table. Each system comes complete with two end frames, two adjustable upper shelves, and one adjustable under table shelf, and a power distribution outlet.

The 12-outlet power strip is available in several power supply outlet types. Shelves fitted with a NEMA-type outlet (PTA261 and PTA271) are shipped from stock. Shelves fitted with SCHUKO (PTA263 and PTA273), or UK (PTA264 and PTA274) type outlets are built upon receipt of order.

Features

- Stable, Free-Standing Design Does not Interfere with Vibration Isolation
- Shelf Height is Adjustable from 43.21" to 70.79" (1.098 to 1.798 m)
- Power Strips Included (Choose from NEMA, SCHUKO, or UK)
- Both Overhead and Under Table Shelves Included
- Shelves Made from Particle Board
- Frame Height: 72.64" (1.845 m)
- Shelf Height (Max): 70.79" (1.798 m)
- Load Capacity (Max): 154 lbs (70 kg) per Shelf
- Shelf Width: 15.75" (400 mm)
- Frame Dimensions: 84.65" x 35.83" (2.15 m x 0.91 m)




Under
Table
Shelf
Included



Power Strip Included.
Four Options Available



SCHUKO
Connector



NEMA
Connector



UK
Connector

Additional Size Options Available

The Free-Standing Shelves featured below are stocked for immediate shipment. Shelves fitted with SCHUKO (PTA263 and PTA273), or UK (PTA264 and PTA274) type outlets are built upon receipt of order. Please visit www.thorlabs.com to view our entire product offering.

Free-Standing Shelf

ITEM #	\$	£	€	RMB	DESCRIPTION
PTA261	\$ 1,385.00	£ 997.20	€ 1,204.95	¥ 11,038.45	Overhead Shelf, 2 x 0.4 m (NEMA)
PTA271	\$ 1,465.00	£ 1,054.80	€ 1,274.55	¥ 11,676.05	Overhead Shelf, 2.5 x 0.4 m (NEMA)